

Why pay more for the same apples?

It's fruitless comparing apples with oranges, but there's nothing wrong with comparing apples with apples, it just makes sense. [Let's see how AlphaApples* compare to Blu-Apples**](#)

* Danville AlphaBite Fast Setting
 ** Parkell Blu-Mousse Super Fast



AlphaApples* vs Blu-Apples**

	Work Time	Set Time	Durometer		Work Time	Set Time	Durometer
	30 Seconds	60 Seconds	Hardness 82		30 Seconds	30 Seconds ¹	Hardness 85 ¹



Other than the variance in set times AlphaBite Fast Setting has been formulated to the same specifications of Blue-Mousse Superfast and has the same handling properties giving you highly stable and accurate results time after time.

There's no comparison when it comes to price!

DAN93352
Alphabite Fast Setting
 2x 50ml Automix Cartridges
 Includes 6 mixing tips



All prices exclude GST

- Also Available:**
- 4 Pack - 4 x 50ml Cartridges
 - Econo Kit - 24 x 50ml Cartridges
 - Super Econo Kit - 96 x 50ml Cartridges



¹ Durometer and set time information taken from Blu-Mousse & Green-Mousse Instructions for Use. www.parkell.com. March 2017

Make the smart switch to AlphaBite and start saving today!