

engler Son-Mate II

Ultrasonic Scaler / Polisher Combination Unit



The SON-MATE 2 is one of the most reliable dental combination units on the market. It is a potent and powerful tool against periodontal disease. At the touch of the switch, you alternate from scaling to polishing plus a variety of other operations.

There are 9 different tips available for the SON-MATE 2. The heat-treated surgical stainless steel tips glide smoothly across the tooth surface while the water flows through the tip DIRECTLY to the operative site, gently cooling and lubricating the tooth whilst washing away the dislodged debris.

- Scaler Hand piece
- 4 Assorted tips
- Water Line with filter and quick disconnect
- 35,000 rpm variable speed precision balanced micro motor
- Forward/reverse switch
- Autoclavable prophylaxis angle for use with snap-on rubber cups
- Autoclavable straight hand piece
- Rubber rest for micromotor
- Transducer (stack) for Scaler Hand piece
- Tip Wrench
- Foot Pedal for Scaler & Polisher
- Digital Power Display
- MicroController based digital auto-tuned circuit
- Quick connect all accessories
- MicroMotor Speed control

ENGA57000
\$3995.00
Usually \$4051.75 / Each

Hardened Glass Lens

90° Direct Light

Stainless Steel Head

360° Rotation

ON/OFF Button

Display Screen

Time Button

Mode Indicators

Mode Button

Charging Dock

premium plus™


DUAL WAVELENGTH LED CURING LIGHT

Light Head **10mm**
Wavelength **390-430nm (lower)**
440-480nm (upper)

Tack	2,000 mW/cm ² For 1 Second±10%
Turbo	2,500 mW/cm ² For 3 Seconds±10%
3'+ 3'	2,000 mW/cm ² For 3 + 3 (1 sec blink) Seconds±10%
6'+ 6'	1,500 mW/cm ² For 6 + 6 (1 sec blink) Seconds±10%
Full	1,200 mW/cm ² For 5, 10, 15, 20 Seconds±10%
Ramp	First 5 sec 0-1,200 mW/cm ² Next 5, 10, 15, 20 sec 1,200 mW/cm ² ±10%
Soft Cure	First 5 sec 0-900 mW/cm ² Next 5, 10, 15, 20 sec 900 mW/cm ² ±10%

Built-in Light Meter

PPIC02-D
\$553.25
Each



DISPOSABLE PROPHY ANGLES

Premium Plus

- Fits all standard straight handpieces.
- Inclined angle 105 degree.
- Guaranteed no dislodging of cup.

Take **20% OFF**

PPI3399-105R/S
\$57.68
Usually \$72.10 / 100 Pack

LATEX FREE PROPHY CUPS

Premium Plus

Take **20% OFF**

Latch / RA PPI222
\$28.84
Usually \$36.05 / 100 Pack

Screw In PPI221
\$24.72
Usually \$30.90 / 100 Pack

Snap On PPI220
\$13.16
Usually \$16.45 / 100 Pack

PERSONALISED RECALL CARDS



THE BUNCH



FREE personalisation on all recall cards!

Pre-printed with your practice details and logo. Templates available or design your own!

100-400 Cards	500-900 Cards	1000+ Cards
SPRPC-04538	SPRPC-04538	SPRPC-04538
\$60.75	\$57.00	\$52.00
/ 100 Pack	/ 100 Pack	/ 100 Pack

CONTRA ANGLE HANDPIECE ATTACHMENT 1:1

Premium Plus



Take
20%
OFF

PPIL20-1L

\$123.60

Usually \$154.50 / Each

STRAIGHT HANDPIECE ATTACHMENT 1:1

Premium Plus



Take
20%
OFF

PPIL10-1

\$123.60

Usually \$154.50 / Each

E-TYPE SNAP ON PROPHY CONTRA ANGLE

Nakamura



Universal E-type contra angle attachment for prophy snap-on.

Features

- Gear ratio 1:1
- Snap-on prophy
- Standard type head
- Fits air or electric driven E-type motor
- Made in Japan

Only
11
LEFT

NAKEC-20PK

\$120.61

Usually \$172.30 / Each

Take
30%
OFF

Have you discovered the Bactol difference?

A growing number of our customers have made the change to Bactol Alcohol Gel.

Here's what our customers are saying about Bactol...

"We've tried numerous others. This is the first gel that hasn't dried my hands out and feels nourishing"

Hendriks Dental

"Nicer on the skin, not as drying & smells better"

Darnell -Manawatu Dental

"The girls in piercing love it!"

Cromwell Pharmacy

"Thank you for allowing us to test this product, our patients prefer Bactol"

Courtney -Feet First Podiatry



Made With Natural
Macadamia Oil

What is the Bactol difference?

- Bactol Won't dry your skin out!
- You do not need to use a moisturiser with it.
- Unlike most other alcohol rubs Bactol has a great cosmetic-quality emollient in it that actually improves the skin!
- The more you use it the better your skin feels.
- Once you rub the emollient in and give it about 30 seconds to dry, it does not leave a sticky residue on your skin like lower quality products.



1L Pod
dispensers
available

WCO020066

\$14.95

/ 500ml Pump Bottle

WCO020059

\$30.00

/ 1 Litre Pod

STERILISATION POUCHES

Healthy Smiles

- 200 per box
- Non-toxic indicators
- Self sealing and reinforced



Everyday
Low
Price!

Size mm	Size Inches	Stock Code	PRICE
57 x 100mm	2 ¼" x 4"	BHMSP57x100	\$5.20
90 x 145mm	3½" x 5¾"	BHMSP90x145	\$8.20
70 x 230mm	2¾" x 9"	BHMSP70x230	\$8.65
90 x 230mm	3½" x 9"	BHMSP90x230	\$13.80
135 x 260mm	5¼" x 10"	BHMSP135x260	\$21.55
190 x 330mm	7½" x 13"	BHMSP190x330	\$34.45
250 x 380mm	10" x 15"	BHMSP250x380	\$77.55

BUY 2 BOXES/ROLLS

OF STERILISATION POUCHES OR STERILISATION ROLLS AND

GET A FREE AUTOCLAVABLE SCRUB BRUSH



STERILISATION ROLL

Healthy Smiles



- Class 3
- Self Sealing
- Steam Sterilisation Only

Size (mm)	Stock Code	Price (Each)
100x200mm	BHMSR100X200	\$36.05
150x200mm	BHMSR150X200	\$48.05



• A true friend reaches for your hand, but touches your heart •

PERSONALISED SYMPATHY CARDS WITH DEEPEST SYMPATHY

Crown Veterinary

- FREE Personalisation your logo + contact details.
- Add your own sympathy message.
- FREE Envelopes.

1-4 10 Packs

CDMSYMP01
\$25.75

5+ 10 Packs

CDMSYMP01
\$22.70

Freephone 0508 333 222

Online www.crownvet.co.nz



New Zealand
Permit No. 1701



Email customerservice@crowndental.co.nz **DECEMBER/JANUARY 2019/2020 EDITION** VALID 01/12/19 – 31/01/20

Fax 03 338 3899

PO Box 640, Christchurch 8140, New Zealand

All prices exclusive of GST and Freight, Errors and omissions excepted (E&OE), Crown Dental + Medical assumes no responsibility or liability for any errors or omissions in the content of this mailer.