

# Opal<sup>®</sup>

## Instrument Grade, High Level Disinfectant



Also Available: **Opal Test Strips** 3 Second Dip Time, Results in 90 secs\*\*

Reuse period up to twice the length of competitors\*

Opal<sup>®</sup> is an Ortho-Phthalaldehyde (OPA) based high level disinfectant that features a reuse period up to twice the length of competitors.\*

Opal<sup>®</sup> is suitable for manual disinfection of compatible semi-critical medical and dental devices and automated reprocessing of Endoscopes and accessories.

With the same contact time and temperature for manual disinfection and automated reprocessing.

- 6-minute contact time
- Used at room temperature (20°C)
- Effective against a wide array of micro-organisms
- Up to 28 day use-life
- Available for manual disinfection
- Ready-to-use, no mixing or activation required
- Contains: Ortho-Phthalaldehyde 5.7g/L
- Minimum Recommended Concentration (MRC) is 0.35% OPA as determined by OPAL Test Strips
- Shelf life is 24 months

\*When the product is decanted from the bottle. \*\*Sold Separately to Opal



WCO150103

**\$128.85**

Usually \$143.16 / 5 Litre

**Whiteley**   
 Medical



## HOSPITAL GRADE DISINFECTANT WIPES

V-Wipes unique formula and high quality apertured fabric allows you to  
**CLEAN & DISINFECT**

- pH neutral – safe for use on most surfaces & equipment.
- Highly absorbent apertured fabric.
- Leaves surfaces streak free.
- Extremely durable.
- Non hazardous.
- Does not contain phenolics or chlorine.
- Non flammable.
- Friendly to use on metals (non corrosive)

WCO210579

**\$6.15**

/ Flatpack of 50

WCO210561

**\$16.40**

/ Cannister of 100

# Try our new LIVE CHAT

Next time you're online at crowndental.co.nz, click on the LIVE CHAT symbol (bottom right of your screen) for prompt answers to you questions.

**Our friendly Customer Service team can answer any questions you have on:**

- Products
- Orders
- Accounts
- Shipping



## UNIVERSAL ABSORBENT SHEET

Tork



- Absorbent sheet helps maintain good hygiene and comfort
- 43cm x 60cm

Everyday  
Low  
Price!

PMO2170057  
**\$52.95**  
/ 100 Pack

## PET HAIR PICK UP ROLLER

SmartPractice



- 60 easy peel sheets
- Removes ticks from clothes
- 29.8 feet x 4 inches per roll
- No backorders

Only  
10  
LEFT

Take  
50%  
OFF

SPR24-607  
**\$6.35**  
Usually \$12.70 / Each



## CORAL PROPHY PASTE

Young

Coral™ Propy Paste is efficient, splatter-free and leaves a high-luster finish.


Available with or without fluoride; spearmint or bubble gum flavors; fine or coarse grit.

Fluoride Option Contains 1.23% fluoride and 0.1 molar phosphate.

Non Fluoride YOU033###  
**\$39.50**  
/ 250gm Bottle

Everyday  
Low  
Price!

Fluoride YOU033###  
**\$42.85**  
/ 250gm Bottle



## ZOOTOPIA STICKERS

Smart Practice

- Set of 6 stickers

SPR292794  
**\$19.80**  
/ 100 Pack

## PLAQCLNZ PET TOOTH-BRUSHES



SmartPractice



- Toothbrushes designed specifically for pets teeth to maximize contact. Tuft count 34 per brush
- Filament diameter .007

SPR24843  
**\$66.65**  
 / 48 Pack

## AUTOCLAVABLE SCRUB BRUSH



Keyport



KMEGIMA25715  
**\$5.40**  
 Each

## LATEX FREE PROPHY CUPS

Premium Plus



Latch / RA  
 PPIIM222  
**\$28.84**  
 Usually \$36.05 / 100 Pack

Screw In  
 PPII221  
**\$24.72**  
 Usually \$30.90 / 100 Pack

Snap On  
 PPII220  
**\$13.16**  
 Usually \$16.45 / 100 Pack

## DISPOSABLE PROPHY ANGLES

Premium Plus

- Fits all standard straight handpieces.
- Inclined angle 105 degree.
- Guaranteed no dislodging of cup.



PPI3399-105R/S  
**\$57.68**  
 Usually \$72.10 / 100 Pack



## Correct Face Mask Selection

On the surface, face masks may look the same but there's a lot more to face masks than meets the eye...

### ASTM Level 1 Low Fluid Resistance

#### SUITABLE FOR

- Patient Exams
- Operatory cleaning / maintenance
- Impressions
- Lab trimming, finishing and polishing
- Orthodontics

#### L1 EAR LOOP FACE MASKS




BHM101-1-1B  
**\$7.75**  
 /Box of 50

### ASTM Level 2 Moderate Fluid Resistance

#### SUITABLE FOR

- Restorative / Composites
- Endodontics
- Prophylaxis
- Sealants
- Scaling or root planning
- Limited oral surgery

#### L2 TIE ON FACE MASKS




BHM143-2-1  
**\$7.75**  
 /Box of 50

### ASTM Level 3 High Fluid Resistance

#### SUITABLE FOR

- Crown preparation
- Implant placement
- Ultrasonic Scaler Use
- Use of Piezo Scaler with water or medicaments
- Periodontal surgery
- Complex oral surgery

#### L3 TIE ON FACE MASKS



BSN36361280  
**\$27.80**  
 /Box of 40

Ask us about MASSIVE bulk discounts **Phone 0508 333 222**

## VISTA-TEC ULTRALIGHT FACE SHIELD KIT

Polydentia



The protective shields are made of transparent, environment-friendly plastic, are scratch resistant and do not yellow when subjected to daylight. They can be used several times as they can be cleaned easily with Polydentia Clean.

Kit Contains 1 Frame + 10 Faceshields  
Available in white frames only.



Take  
**15%**  
OFF

POL5624/W

**\$116.40**

Usually \$136.95 / Per Kit

## STERILISATION POUCHES

Healthy Smiles

- 200 per box
- Non-toxic indicators
- Self sealing and reinforced



Everyday  
Low  
Price!

Size mm	Size Inches	Stock Code	PRICE
57 x 100mm	2 ¼" x 4"	BHMSP57x100	<b>\$5.20</b>
90 x 145mm	3½" x 5¾"	BHMSP90X145	<b>\$8.20</b>
70 x 230mm	2¾" x 9"	BHMSP70x230	<b>\$8.65</b>
90 x 230mm	3½" x 9"	BHMSP90x230	<b>\$13.80</b>
135 x 260mm	5¼" x 10"	BHMSP135x260	<b>\$21.55</b>
190 x 330mm	7½" x 13"	BHMSP190x330	<b>\$34.45</b>
250 x 380mm	10" x 15"	BHMSP250x380	<b>\$77.55</b>

*Good friends are impossible to forget*



## PERSONALISED SYMPATHY CARDS HORSES AT SUNSET

Crown Veterinary

- FREE Personalisation your logo + contact details.
- Add your own sympathy message.
- FREE Envelopes.

View the full range  
at [crownvet.co.nz](http://crownvet.co.nz)

1-4 10 Packs

CDMSYMP08

**\$25.75**

5+ 10 Packs

CDMSYMP08

**\$22.66**

Freephone 0508 333 222

Online [www.crownvet.co.nz](http://www.crownvet.co.nz)



New Zealand  
Permit No. 1701



Email [customerservice@crowndental.co.nz](mailto:customerservice@crowndental.co.nz)

Fax 03 338 3899

SEPTEMBER 2019 EDITION VALID 01/09/19 – 30/09/19

PO Box 640, Christchurch 8140, New Zealand

All prices exclusive of GST and Freight, Errors and omissions excepted (E&OE), Crown Dental + Medical assumes no responsibility or liability for any errors or omissions in the content of this mailer.